

5641. Single cell spatial molecular imaging of 119-plex proteins in clinical cancer samples in response to personalized treatment

Zachary R. Lewis¹, Giang Ong¹, Tièn Phan-Everson¹, Evelyn Metzger¹, Brian Birditt¹, Emily Brown¹, Kan Chantranuvatana¹, Chris Corless², Brian Filanoski¹, Gary Geiss¹, Evie Hobbs², Brett Johnson², Taylor Kelley², Mithra Korukonda¹, Charles Lopez², Rhonda Meredith¹, Anastasiya Olson², Erin Piazza², Jason Reeves¹, Alyssa Rosenbloom¹, Kiara Siex², Hye Son Yi¹, Edward Zhao¹, Joseph M. Beechem¹, Gordon Mills²

1. NanoString® Technologies, Seattle WA. 2. OHSU, Knight Cancer Institute, Portland OR

Abstract

The power of spatial biology lies in the integration of multiple scales of information from subcellular to tissue scale. Until recently, spatial analysis of protein biomarkers in tissues was limited to a few markers at a time using traditional IHC colorimetric or fluorescent readout. Here, we demonstrate the ability of NanoString's CosMx™ Spatial Molecular Imager (SMI) platform to quantify more than 100 proteins encompassing key targets in immuno-oncology and tumor biology, localize the proteins and analyze protein expression at a single cell level. Key to the technology is the use of fully automated fluidics and imaging systems, short turnaround time, and high sensitivity. The CosMx protein assay has been optimized for FFPE samples, which represent the largest collection of biospecimens available for clinical investigation. To that end, we explored the utility of CosMx spatial proteomics on a series of clinical samples from cancer patients, including the Serial Measurements of Molecular and Architectural Response to Therapy (SMMART) Clinical Trials Program, across multiple cancer types.

CosMx SMI's high-level protein multiplexing capabilities enabled spatial analysis of metastatic tumors in response to personalized treatment for a single patient over time. Detection of phosphoproteins also allows for analysis of the impact of kinase inhibitor treatment on the spatial environment in longitudinal biopsies. The highly-multiplexed spatial analysis of proteins in longitudinal metastatic breast cancer biopsies under therapeutic pressure provides a unique opportunity to understand evolution of tumors and develop and implement therapeutic approaches that can directly target mutations arising in the tumor cells while effectively engaging the immune system.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

Highlights

- Spatial profiling of **14 breast cancer patients** undergoing personalized treatment
- **33 biopsies** (FFPE core needle) across 19 slides
- **119-plex** CosMx protein assay
- **Post-translational** modifications captured
- Immune **cell composition and expression** characterized over treatment

Longitudinal study of patients undergoing personalized treatment for breast cancer

Serial Measurements of Molecular and Architectural Responses to Therapy (SMMART®) Clinical Trials Program

The patients are followed longitudinally, which in this case means that they are monitored before, during and after being given study treatments.

The treatments are tailored to each individual patient in the study according to the characteristics of the individual's cancer and genetics.

Adaptive treatment in real time. Therapy for each individual patient is tailored, then re-evaluated, and changed during the trial according to changes observed in the patient's response.

Objectives:

- Identify new treatments for cancer that last longer (are more **durable**)
- Allow better quality of life (are more **tolerable**) for patients with advanced disease
- Understand why therapies often stop working

CosMx Spatial Molecular Imaging

CosMx protein assays detect antibody localization by hybridization of a bright, photocleavable reporter construct. The reporter construct can be gently removed by photocleavage. The commercial 64-plex Human Immuno-oncology panel was combined with conjugated antibodies for an additional 51 targets, as well as a four-target segmentation cocktail, yielding a 119-plex panel.

