

Rudy Van Eijsden¹, Alyssa Rosenbloom¹, Shilah Bonnett¹, Christine Kang¹, Mark Conner¹, Megan Vandenberg¹, Erin Piazza¹, Brian Filanoski¹, Rhonda Meredith¹, Hye Son Yi¹, Lori Hamanishi¹, Eduardo Ignacio¹, Nadine Nelson², Michael Prater², Vik Devgan¹, Melanie Moon¹, Lesley Isgur¹, Terence Theisen¹, Margaret Hoang¹, Gary Geiss¹, and Joseph M. Beechem^{1*}

¹NanoString® Technologies, Seattle, WA, USA, ²abcam plc, Discovery Drive, Cambridge, UK


Abstract

The advancement of spatially resolved, multiplex proteomic and transcriptomic technologies has revolutionized and redefined the approaches to complex biological questions pertaining to tissue heterogeneity, tumor microenvironments, cellular interactions, and therapeutic response. While spatial transcriptomics has traditionally led the way in plex, multiple studies have demonstrated a poor correlation between RNA expression and protein abundance, owing to transcriptional and translational regulation, target turnover, and most critically, post-translational protein modifications. Therefore, a more holistic, ultra-high-plex, and high-throughput proteomic atlas approach becomes critical for the next phase of discovery biology. Here, we present a barrier-breaking, spatial proteomics panel that was designed to accelerate scientific discoveries.


A Digital Spatial Profiler platform is uniquely suited to support high-plex proteomics, allowing for the simultaneous analyses of proteins from discrete regions of interest (ROIs) in FFPE tissue sections while preserving spatial context. The assay relies upon abcam antibodies coupled to photocleavable DNA barcodes readout with NGS sequencing, allowing for theoretically unlimited plex. Here we introduce the Human Immuno-Oncology Proteome Atlas (IPA), a 570+ antibody-based proteomic discovery panel, compatible with immunohistochemistry on FFPE tissues with NGS readout. IPA is the highest-plex multi-omic (~610-plex proteins and >18,042 genes) study ever implemented for spatial biology. When we compared the diseased tissue to normal tissue, we observed an upregulation of specific pathways associated with tumorigenesis and inflammation. Furthermore, we observed distinct differences in proteomic and transcriptomic landscape between pathologies. The cutting-edge, data-driven, expert-curated IPA panel is at the forefront of spatial proteomics, empowering the researcher for the acceleration of biological discoveries.

Here we demonstrate the performance of IPA on various cell lines and tissue. Additionally, we show the power of IPA, using the spatial multi-omic assay along with the GeoMx® Whole Transcriptome Atlas (~18,000 transcripts), a 30-plex custom antibody panel and microbiome-curated RNA custom spike-in (~42 transcripts) to evaluate 70 different colon disease samples across 4 pathologies including adenocarcinoma, hyperplasia, and chronic inflammation. This is the highest-plex multi-omic (~610-plex proteins and >18,042 genes) study ever implemented for spatial biology. When we compared the diseased tissue to normal tissue, we observed an upregulation of specific pathways associated with tumorigenesis and inflammation. Furthermore, we observed distinct differences in proteomic and transcriptomic landscape between pathologies. The cutting-edge, data-driven, expert-curated IPA panel is at the forefront of spatial proteomics, empowering the researcher for the acceleration of biological discoveries.

Spatial Proteogenomics reveals reduced phosphorylation of Cytokeratin 8 in Malignant Adenocarcinoma compared to chronically Inflamed colon




GeoMx® Digital Spatial Profiler Workflow



GeoMx Immuno Oncology Proteome Atlas is fully compatible with existing and well validated GeoMx Digital Spatial Profiler workflows


Spatial Proteogenomic Workflow #2: Seq Code RNA Assay and Pro Code Protein Assay



Designing an Immuno Oncology Proteome Atlas

Curated and validated IO Content

- Use Clinical Proteomic Tumor Analysis Consortium (CPTAC) data
- Focus on post-translational modifications
- Use high-quality abcam antibodies
- Employ Immuno Oncology subject matter experts
- Human Specific
- 570+ Proteins, Mapped to 556 Unique Genes
- 77 Functional Annotations
- All Hallmarks of Cancer



GeoMx Immuno Oncology Proteome Atlas Functional Pathways

GeoMx Immuno Oncology Proteome Atlas Public Dataset

Visit
GeoMx® IO Proteome Atlas | NanoString

Spatial Proteogenomics on 7 and 8 week Human Embryo shows distinct temporal organ differentiation




Figure 5: GeoMx DSP images of 7 wk (a) & 8 wk (b) Hu embryo (*donated tissue from ectopic pregnancy). c) UMARS of gene expression (RNA) & protein expression in 7 and 8 wk embryonic organ systems. SpatialOmics Overlay comparisons of key Protein & RNA targets in 7 wk (d) & 8 wk (e) embryo slices. f) Box Plots of key protein & RNA target expression across 7 and 8 wk organ systems.

*BS de Bakker, MJA van den Hoff et al. Dutch Fetal Biobank, Amsterdam UMC.

References

- Merritt CR, et al. Multiplex digital spatial profiling of proteins and RNA in fixed tissue. *Nat Biotechnol*. 2020;38(5):586-599.
- Zollerger DR, et al. GeoMx® RNA Assay: High Multiplex, Digital, Spatial Analysis of RNA in FFPE Tissue. *Methods Mol Biol*. 2020;2148:331-345.
- Zimmerman SM, et al. Spatially resolved whole transcriptome profiling in human and mouse tissue using Digital Spatial Profiling [published online ahead of print, 2022 Sep 13]. *Genome Res*. 2022;gr.276263.11.
- Bonnett S, A., Rosenbloom A, et al. Ultra High-Plex Spatial Proteogenomic Investigation of Glioblastoma Multiforme Immune Infiltrates Reveals Distinct Protein and RNA Expression Profiles. (2023) *Cancer Research Communications*. 3(5):763-779. <https://doi.org/10.1158/2767-9764.CRC-22-0396>
- <https://nanostring.com/products/geo-mx-digital-spatial-profiler/io-proteome-atlas/>
- www.abcam.com

Scan here to download or learn more

