


Alyssa Rosenbloom¹, Ashley Heck¹, Kimberly Young¹, Aster Wardhani¹, Max Walter¹, Giang Ong¹, Rustem Khafizov¹, Edward Zhao¹, Sangsoo Woo¹, Lidan Wu¹, Mithra Korukonda¹, Zachary Lewis¹, Tien Phan-Everson¹, Shilah A. Bonnett¹, Claire Williams¹, Patrick Danaher¹, Rachel Liu¹, Jose Ceja Navarro¹, Christine Kang¹, Melanie Moon¹, Jim Forrester¹, Trieu My Van¹, Jason Reeves¹, Gary Geiss¹, David W. Ruff¹, Joseph M. Beechem¹

Introduction

68 plex human protein panel coverage


The Spatial Molecular Imager (SMI) addresses this problem by enabling highplex detection of both proteins and RNAs on the same slide to improve brain cell segmentation and assignment of transcripts to single cells. In this multi-omic system, 68 proteins and over 6,000 transcripts covered by the Human 6K Discovery Panel are targeted on the exact same FFPE human brain section



Multi-omic workflow overview




Sample preparation, flow cell assembly and target detection

Decoded RNA targets + protein
cell segmentation

Cell segmentation based on protein run morphology captures neuron and glial processes in addition to cell bodies and associated transcripts



“New” cell segmentation framework output



Protein data-based unsupervised cell typing Protein markers RNA markers Protein data-based cell typing +



Protein data-based unsupervised cell typing		Protein markers	RNA markers	Unsupervised neighborhood analysis identifies tissue
Spatial	UMAP			



By combining high-plex neuro protein and RNA assays, we gain extraordinary access to the spatial dynamics of multimodal interactions. Specifically, the high-plex SMI multi-omic assay:

- Covers 68 proteins focused on neural cell typing and neural pathology/ inflammation.
- Profiles >4,900 neuroscience-related genes that cover areas ranging from basic neurobiological signaling mechanisms to neurodegenerative disease pathology on the exact same slide.
- Enables unparalleled segmentation of extended cell processes for neurons and glia.
- Rescues many of the 50% unassigned orphan RNA transcripts and quantitates the sub-cellular spatially localized expression in the human brain, an area of extensive biological activity.
- Enables refined cell typing that takes advantage of both protein and RNA data.
- Opens numerous avenues for biological inquiry by combining cell typing with neighborhood analyses.



The CoSeq™ SMI and decoder probes are not offered and/or delivered to the following UPC member states* for use in these countries for the detection of RNA in a method used for the detection of a plurality of analytes in a cell or tissue sample without the consent of the President and Fellows of Harvard College (Harvard Corporation) as owner of the Unitary Patent EP 1 08 782 B1. The use for the detection of RNA is prohibited without the consent of the President and Fellows of Harvard College (Harvard Corporation).

*Austria, Belgium, Bulgaria, Denmark, Estonia, Finland, France, Germany, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Portugal, Slovenia, Sweden.

The CoSeq™ SMI and decoder probes are not offered and/or delivered to the Federal Republic of Germany for use in the Federal Republic of Germany for the detection of cellular RNA, messenger RNA, microRNA, ribosomal RNA and any combinations thereof in a method used in fluorescence in situ hybridization for detecting a plurality of analytes in a sample without the consent of the President and Fellows of Harvard College (Harvard Corporation) as owner of the Unitary Patent EP 1 08 782 B1. The use for the detection of cellular RNA, messenger RNA, microRNA, ribosomal RNA and any combinations thereof is prohibited without the consent of the President and Fellows of Harvard College (Harvard Corporation).

The CosMax™ SMI and decoder probes are not offered and/or delivered to the Federal Republic of Germany for use in the Federal Republic of Germany for the detection of cellular RNA, messenger RNA, microRNA, ribosomal RNA and any combinations thereof in a method used in fluorescence in situ hybridization for detecting a plurality of analytes in a sample without the consent of the President and Fellows of Harvard College (Harvard Corporation) as owner of the German part of EP 2 794 928 B1. The use for the detection of cellular RNA, messenger RNA, microRNA, ribosomal RNA and any combinations thereof is prohibited without the consent of the use of the President and Fellows of Harvard College (Harvard Corporation).